## **Tracheostomy Tube**

### **Model:**

SuctionPlus Tracheostomy Tube with Evacuation Lumen

### **DESCRIPTION**

Tracheostomy Tube is a device used in general anesthesia, intensive care and emergency medicine for airway management and mechanical ventilation. The tube is inserted into a patient's trachea through an incision in the trachea (the windpipe).

SuctionPlus Tracheostomy Tube with Evacuation Lumen is applicable for child, adult patients.

SuctionPlus Tracheostomy Tube with Evacuation Lumen consists of main tube, neck-plate, cuff, connector, suction tube, suction connector, inflation tube, pilot balloon, non-return valve, introducer, ribbon, heat shrinkable film and cover.

### **FEATURES**

- Full length radio-opaque line for detection of correct placement by X-ray.
- ISO standard connector for universal connection to ventilation equipment
- Printed neck plate with size information for easy identification.
- Ribbon provided in the pack for fixation of the tube.
- Smooth rounded tip of the obturator reduced trauma during insertion.
- Rigid blister pack provides maximum protection to the tube.

### **INTENDED USE**

The Tracheostomy Tube is used in general anesthesia, intensive care and emergency medicine for airway management and mechanical ventilation. The tube is inserted into a patient's trachea through an incision in the trachea (the windpipe).

### **INDICATIONS**

The indications are expanded with the development and widespread application of tracheotomy, mainly to keep the respiratory tract unimpeded for the patients unable to receive oropharyngeal intubation or nasotracheal intubation, specially including: ① upper respiratory tract obstruction: upper respiratory tract obstruction caused by foreign body in respiratory tract, diphtheria, acute laryngopharyngitis and other infections, trauma, angioneurotic edema, oropharyngeal neoplasm, tracheomalacia, neck surgery or radiation-induced scar contracture, obstructive sleep apnea syndrome or any other reasons; ② Retention of secretions in the lower respiratory tract: retention of secretions in the lower respiratory tract caused by craniocerebral trauma, intracranial or peripheral neuropathy, disturbance of consciousness, tetanus, inhalation burn, post-major thoracic and abdominal surgery, bronchiectasis, chronic pneumonia, laryngeal dysfunction, hypofunction or loss of cough and sputum excretion, and aggregation of secretions; ③ Central or peripheral respiratory depression: central or peripheral respiratory depression caused by encephalitis, myelitis and other infections, poisoning, high fever, high intracranial pressure, cerebral hernia, traumatic brain and spinal cord injury, drug inhibition, bilateral vocalcord paralysis, Guillain-Barre Syndrome, amyotrophic lateral sclerosis and myasthenic crisis; ④ Cardiopulmonary insufficiency: patients with severe pulmonary heart disease and poor pulmonary function; ⑤ Prolonged mechanical ventilation: providing stable airway and avoiding laryngeal injury to facilitate the intratracheal administration and to improve the patient comfort; ⑥ Preventive tracheotomy: preventive tracheotomy for keep the respiratory tract unimpeded after surgery of severe trauma of head, face and neck, edema and inhalation injury, oral and maxillofacial surgery and major surgery of ENT, head and neck.

### **CLINICAL BENEFITS**

By inserting into a patient's trachea through an incision in the trachea (the windpipe) to mainly keep the respiratory tract unimpeded for the patients unable to receive oropharyngeal intubation or nasotracheal intubation, it is used in general anesthesia, intensive care and emergency medicine for airway management and mechanical ventilation. The clinical benefit associated with the device is define as direct, patients who can benefit from the device includes: Child, adult patients.

### **PATIENT TARGET GROUP**

SuctionPlus Tracheostomy Tube with Evacuation Lumen: Child, adult patients.

### **INTENDED USER**

Must be operated by trained professionals

### **CONTRAINDICATIONS**

Care must be taken to avoid contact of a LASER beam or an electrosurgical active electrode with this and other tracheostomy tubes. Such contact, especially in the presence of oxygen-enriched mixtures, could result in rapid combustion of the tube with harmful thermal effects and with the emission of corrosive and toxic combustion products, including hydrochloric acid (HCL).

### **ADVERSE REACTIONS**

Standard textbooks and scientific literature should be consulted for any specific adverse reaction information.

### **DIRECTIONS FOR USE**

- Take out the tracheostomy tube from the package.
- Test the cuff by inflation using a Lure-tip syringe.
- Completely deflate cuff after test inflation.
- With the Introducer in place (if supplied), insert the tracheostomy tube through the stoma and position in the tracheal in accordance with currently accepted medical technique.
- Following insertion, remove the Introducer (if used) and slowly inflate the cuff with the minimum amount of air required for effective seal.
- Remove the syringe from the non-return valve.
- Attach the swivel adapter to the 15 mm termination of the tracheostomy tube and connect the ventilator circuit to the swivel adapter. The anti-disconnect lugs on the swivel should be utilized to help against accidental disconnect.
- Secure the tracheostomy tube.
- If using Suction Plus model, can connect the suction port to syringe without needle or suction instruments to proceed suction of secretion.

### **MRI Safety Information**

The Tracheostomy Tube has not been evaluated for safety in the MR environment. It has not been tested for heating or unwanted movement in the MR environment. The safety of Tracheostomy Tube in the MR environment is unknown. Performing an MR exam on a person who has this medical device may result in injury or device malfunction.

### **WARNINGS**

- Read all warnings and instructions before use. Improper use can result in serious or fatal illness or injury.
- This device can only be used by trained professionals.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Do not use if it's damaged or irregularly shaped.

- Do not use after the expiry date.
- For single use only.

## PRECAUTIONS

- If cuffed Tracheostomy Tube is used, test the cuff to ensure no leakage and inflate with minimum necessary volume of air. And retain the cuff pressure between 25 to 30 cmH<sub>2</sub>O and monitor the cuff pressure 3 times a day.
- The product cannot retain in patient for more than 29 days. Medical care person need to do nursing of the patient's tracheostomy and judge when to do the extubation according to the patient's actual situation.
- Inflation of cuff by "feel" only or by using a measured amount of air is not recommended since resistance is an unreliable guide during inflation.
- Deflate cuff prior to repositioning the tube. Movement of the tube with cuff inflated could result in patient injury or damage to the cuff, and will require a tube change.
- Diffusion of nitrous oxide (N<sub>2</sub>O) mixture, oxygen or air from a surrounding gas mixture may either increase or decrease cuff volume and pressure. In order to limit such effect it is recommended to inflate the cuff with the same gas mixture, which will contact its external surface.
- Do not overinflate cuff. Over inflation can result in tracheal damage, rupture of the cuff with subsequent deflation, or in cuff distortion which may lead to airway blockage.
- Do not use needle syringe to inflate cuff. Use Luer syringe.
- In the event, abnormal anatomy is encountered or unusual positioning of the head or neck is required, care should be exercised to avoid the potential for kinking of the tracheostomy tube. The use of a reinforced tube should be considered for unusual cases.
- Attach the connector firmly in both the tracheostomy tube and adapter on the ventilation equipment to prevent disconnection during use.
- Non-standard dimensions of some connectors on ventilators or anesthesia equipment may affect secure mating with the Tracheostomy Tube 15mm connector. Use only with equipment featuring standard 15mm connectors.
- Expert clinical judgment should be executed in the selection of the appropriate size tracheostomy tube for each patient.
- Intubation and extubation should be performed following currently accepted medical techniques.
- Each tube's cuff, pilot balloon, and valve should be tested by inflation prior to use.
- Care should be taken to avoid damaging the thin-walled cuffed during intubation.
- Avoid using Lidocaine Topical Aerosol as it has been associated with the formation of pinholes in PVC cuffs (Jayasuriya, K.D., and Watson, W.F.: PVC. Cuffs and Lignocaine-based Aerosol. Brit. J. Ann. 53:1368, 1981). The same authors report that the lidocaine hydrochloride solution does not have this effect.
- Tree-way stopcocks or other devices should not be left inserted in the inflation valve for extended periods of time. The resulting stress may break the valve housing and cause the cuff to deflate.
- Alcohol-based solutions should not be in contact with the tube.
- Excessive use of lubricant may cause blurring or blockage of the tube.
- Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.

- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the user's place.
- The maximum withstand negative pressure of the suction tube is 40kpa.

### STORAGE CONDITIONS

- Store in a cool and dry place, keep away from sunlight.
- Protect product from moisture and excessive heat.
- Avoid prolonged exposure to ultraviolet, sunlight and fluorescent light.
- Store in manner preventing crushing.

### The Key Performance Characteristics of the Tracheostomy

#### Key Performance 1: Connector

The tracheostomy tube could be fitted, at the machine end, with a conical connector having a 15 mm cone complying with ISO 5356-1.

#### Key Performance 2: **Cuff Performance** (for cuffed type)

SuctionPlus Tracheostomy Tube with Evacuation Lumen:



Tube size (mm)	Minimum trachea diameter (mm)	Maximum trachea diameter (mm)	Cuff pressure hPa(cmH2O)	Leakage rate range (mL/h)			
				Minimum trachea		Maximum trachea	
				50th percentile	90th percentile	50th percentile	90th percentile
6.0	14.5	19.0	27	270.2	707.4	56.5	207.0
6.5	15.2	19.0	27	413.2	838.9	57.1	307.3
7.0	15.7	20.0	27	490.8	866.5	0.0	1.8
7.5	16.1	20.0	27	501.4	1162.3	0.0	94.4
8.0	16.1	20.0	27	673.3	1418.6	186.0	471.4
8.5	17.0	20.0	27	801.8	1816.4	59.5	774.3
9.0	17.1	20.0	27	666.8	1378.5	552.5	1354.7

#### NOTE:

\* Per ISO 5366.

\*The performance information shown above was collected using a bench test that is intended to provide a comparison of the sealing characteristics of the cuffs only in a laboratory setting. The bench test is not configured or intended to predict performance in the clinical setting.

## MEANING OF SYMBOLS ON PACKAGE

	<b>Date of manufacture</b>		<b>Use-by-date</b>		<b>Batch code</b>
	<b>Do not re-use</b>		<b>Consult instructions for use</b>		<b>Do not use if package is damaged</b>
	<b>Sterilized using ethylene oxide</b>		<b>Do not resterilize</b>		<b>Single sterile barrier system</b>
	<b>Medical Device</b>		<b>Manufacturer</b>		<b>Authorized representative in the European Community</b>
	<b>CE Marked Product</b>		<b>Keep dry</b>		<b>Keep away from sunlight</b>
	<b>Importer</b>		<b>This way up</b>		<b>Fragile, handle with care</b>
	<b>Catalogue number</b>		<b>Not contains or presence of phthalate.</b>		<b>Temperature limits: -15°C~49°C</b>

WELL LEAD MEDICAL CO., LTD.

Address: C-4 Jinhu Industrial Estate, Hualong, 511434 Panyu, Guangzhou, PEOPLE'S REPUBLIC OF CHINA

E-mail: [info@welllead.com.cn](mailto:info@welllead.com.cn) [www.welllead.com.cn](http://www.welllead.com.cn)

Tel: +86 20 84758878

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg, GERMANY

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