Intubating Stylet

DESCRIPTION

The Intubating Stylet is a malleable, aluminum stylet covered with PVC sleeve. The Intubating Stylet provides a positive aid to intubation. It is intended to facilitate quick curving of the Tracheal Tube in preparation for intubation. Allow the Tracheal Tube to be directed more easily, especially in the difficult intubation.

The Intubating Stylet consists of PVC sleeve and malleable aluminium
The product can be used on neonate, infant, child, adult patients.

Product	Type	Size
Intubating Stylet	n type with sleeve	6Fr, 8Fr, 10Fr, 12Fr, 14Fr
	n type without sleeve	6Fr, 10Fr, 14Fr
	e type with sleeve	6Fr, 10Fr, 14Fr

FEATURES

- One piece with smooth end.
- In-built aluminum rod, wrapped with clear PVC.
- Malleable aluminium which allows the stylet to be easily made to any desired shape.
- Latex free.

INTENDED USE

Intubating Stylet is intended to shape appropriate curve of endotracheal tubes prior intubation to facilitate insertion of endotracheal tubes.

INDICATIONS

It is for disposable use by the medical department when auxiliary insertion for tracheal intubation.

CLINICAL BENEFITS

Increasing the tracheal insertion success rate and shortening tracheal insertion time so as to improve patients operative quality.

PATIENT TARGET GROUP

Neonate, infant, child, adult patients.

INTENDED USER

Must be operated by trained professionals

CONTRAINDICATIONS

There are no known contraindications.

ADVERSE REACTIONS

- Accidental extubation
- Aspiration of detached tip or sheared plastic sheath
- Damage to pyriform fossa, cricothyroid membrane, or cricopharyngeal membrane with resulting subcutaneous emphysema, mediastinitis or pneumothorax
- Lacerations of trachea or larynx and tissue trauma

DIRECTIONS FOR USE

- The Intubating Stylet is used with oral Tracheal Tube.
- Select the suitable size of Intubating Stylet.
- Peel off the package, generously lubricate the stylet and insert the stylet into the Endotracheal Tube from the machine end. Make the Tracheal Tube with stylet to desired shape as needed. And then Insert the Endotracheal Tube with Intubating Stylet into trachea, and make Tracheal Tube successfully enter into trachea during intubating.
- After Tracheal Tube is intubated, take out the Intubating Stylet.

WARNINGS

- Read all warnings and instructions before use. Improper use can result in serious or fatal illness or injury.
- This device can only be used by trained professionals.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Do not use if it's damaged or irregularly shaped.
- Do not use after the expiry date.
- For single use only.

PRECAUTIONS

- Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Care should be taken not to abrade the plastic sheath on the edges of the 15 mm connector during the insertion and removal of the stylet.
- This device must not be cut to length by the user.
- Extension of the stylet tip beyond the distal tip of the tracheal tube could result in tissue trauma.
- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the user's place.

STORAGE CONDITIONS

- Store in a cool and dry place, keep away from sunlight.
- Protect product from moisture and excessive heat.
- Avoid prolonged exposure to ultraviolet, sunlight and fluorescent light.
- Store in manner preventing crushing.

The Key Performance Characteristics of Intubating Stylet

Shaping Performance

When inserting the Intubating Stylet into a tracheal tube according to the table 1, 2, 3, the tracheal tube shall be able to maintain its shape after being arbitrarily shaped and placed for 1h.

Table 1 Size of Intubating Stylet (n type with sleeve) inserted into tracheal tube

Size of Intubating Stylet	6Fr	8Fr	10Fr	12Fr	14Fr
Size of applied tracheal tube(ID)	2.5-5.5	3.5-6.5	4.0-7.5	≥5.0	≥5.5


Table 2 Size of Intubating Stylet (n type without sleeve) inserted into tracheal tube

Size of Intubating Stylet	6Fr	10Fr	14Fr
Size of applied tracheal tube(ID)	2.5-5.5	4.0-7.5	≥ 5.5

Table 3 Size of Intubating Stylet (e type with sleeve) inserted into tracheal tube


Size of Intubating Stylet	6Fr	10Fr	14Fr
Size of applied tracheal tube(ID)	3.0-5.5	4.0-7.5	≥ 5.5

MEANING OF SYMBOLS ON PACKAGE

	Date of manufacture		Use-by-date		Batch code
	Do not re-use		Consult instructions for use		Do not use if package is damaged
	Sterilized using ethylene oxide		Do not resterilize		Single sterile barrier system
	Medical Device		Manufacturer		Authorized representative in the European Community
	CE Marked Product		Keep dry		Keep away from sunlight
	Importer		This way up		Fragile, handle with care
	Not contains or presence of phthalate.		Catalogue number		Temperature limits: -15°C~49°C

WELL LEAD MEDICAL CO., LTD.
 Address: C-4 Jinhua Industrial Estate, Hualong, 511434 Panyu, Guangzhou, PEOPLE'S
 REPUBLIC OF CHINA
 E-mail: info@welllead.com.cn www.welllead.com.cn
 Tel: +86 20 84758878

Shanghai International Holding Corp. GmbH (Europe)
 Address: Eiffestrasse 80, 20537 Hamburg, GERMANY

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