Suction Catheter

DESCRIPTION

Suction Catheters are used to remove sputum and secretion from respiratory tract. The suction catheter consists of suction tube and connector. The tube is made from medical grade PVC and the connector is made from PVC or K-Resin (depends on the type of connector).

This product can be used on patient who need to suction of secretion.

FEATURES

- Color-coded connector for easy size identification
- Soft distal end facilitates comfortable insertion
- Different type of connectors available

INTENDED USE

Suction Catheter is used to suck the sputum and secretion in respiratory tract, in order to prevent plugging of the airways. The catheter is used by directly inserted into the throat or through inner lumen of Tracheal Tube.

CLINICAL BENEFITS

The Suction Catheter can suck the sputum and secretion in respiratory tract, in order to prevent plugging of the airways.

The clinical benefit associated with the device is define as direct, it is intended used for neonate, infant, pediatric and adult patients who need aspiration of sputum or secretion.

PATIENT TARGET GROUP

Neonate: 6Fr

Infant:6Fr、 8Fr

Child: 6Fr、 8Fr、 10Fr

Adult:10Fr、 12Fr、 14Fr、 16Fr、 18Fr、 20Fr

INTENDED USER

Must be operated by trained professionals

INDICATIONS

Clinically, it is often used in patients who cannot effectively cough and expel phlegm caused by critical illness, old age, coma and anesthesia.

CONTRAINDICATIONS

- Nasal airway obstruction
- Obvious deviation of nasal septum
- Basilar skull fracture
- Fracture of nasal bone
- Abnormal coagulation mechanism

ADVERSE REACTIONS

- Airway mucosal injury
- Aggravated hypoxia
- Pulmonary atelectasis
- Bronchospasm
- Arrhythmia
- Changes in blood pressure (elevated blood pressure)

DIRECTIONS FOR USE

- Wear gloves, mask and take a face shield.
- Choose the appropriate models and sizes. Check the model and size on the product label to make sure the proper model and size is used.
- Check if the package is in good condition and make sure product is intact.
- Take out the Suction Catheter from the packaging and put in the sterile field.
- Connect the connector of the Suction Catheter to the suction device.
- Turn on the suction, check the suction ability and adjust the vacuum.
- Disconnect the patient from the ventilation source.
- Insert the Suction Catheter into the oral cavity or Tracheal Tube to perform suction. After suctioning, connect the ventilation source again to the Tracheal Tube in use.

PRECAUTIONS

- Expert clinical judgment should be exercised in the selection of the appropriate size tube for each individual patient.
- Do not force tube if any obstruction is encountered during insertion.
- Each suction time should be less than 15 seconds to prevent hypoxia.

WARNINGS

- Read all warnings and instructions before use. Improper use can result in serious or fatal illness or injury.
- This device can only be used by trained professionals.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Do not use if it's damaged or irregularly shaped.
- Do not use after the expiry date.
- For single use only.
- Observe the situation of the patients, stop suction and provide oxygen if there is any significant change on heart rate, blood pressure, respiration or blood oxygen saturation.
- Suck the oral cavity secretion first and then suck the secretion inside the Tracheal Tube. If patient are with Tracheostomy Tube, first suck the opening of Tracheostomy Tube and then suck the secretion in oral and nasal cavity.
- Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the user's place.

- For suction catheters intended to be used with tracheal tubes, the instructions for use shall contain a warning to the effect that the suction catheter is to be fully removed from the tracheal tube before the tracheal tube is cut to length.

STORAGE CONDITIONS

- Store in a dark, cool and dry place.
- Protect product from moisture and excessive heat.
- Avoid prolonged exposure to ultraviolet, sunlight and fluorescent light.
- Store in manner preventing crushing.

APPENDIX A The Key Performance Characteristics of Suction Catheter


- Shaft performance

The shaft shall not collapse when subjected to a pressure ≤ 40 kPa below atmosphere for ≤ 15 s at a temperature of (23 ± 2) °C with the tip occluded and if present, the suction control device set to maximum.

- Suction control device performance

The residual vacuum of the suction catheter, if provided with a suction control device, shall be < 0.33 kPa (3.4 cmH₂O). Check conformance by the tests given in Annex C according to ISO 8836-2019.

MEANING OF SYMBOLS ON PACKAGE

	Date of manufacture		Use-by-date		Batch code
	Do not re-use		Consult instructions for use		Do not use if package is damaged
	Sterilized using ethylene oxide		Do not re-sterilize		Single sterile barrier system
	Medical Device		Manufacturer		Authorized representative in the European Community
	CE Marked Product		Keep dry		Keep away from sunlight
	Importer		This way up		Temperature limits: -15°C~49°C

Fragile, handle
with care

Not contains or
presence of
phthalate.
(if applicable)

Catalogue
number


WELL LEAD MEDICAL CO., LTD.

Address: C-4 Jinhua Industrial Estate, Hualong, 511434 Panyu, Guangzhou, PEOPLE'S
REPUBLIC OF CHINA

E-mail: info@welllead.com.cn www.welllead.com.cn

Tel: +86 20 84758878

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg, GERMANY

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