Urethral Dilator

Model

Hydrophilic Coating Dilator

Hydrophilic Coating S Shape Dilator

DESCRIPTION

The Urethral Dilator is used for urethral dilation. The Urethral Dilator has hydrophilic coating, and is offered in French size ranging from 8Fr to 26Fr. It can be provided in set form (includes any size of dilator) or single dilator.

It is supplied sterile and intended for single use only.

Stainless steel guide wire, PTFE coating guide wire or PTFE guide wire are optional.

The product can be used on urethral strictures patient.

FEATURES

- Made of polyethylene
- Smooth surface for less trauma and easy advancement
- Tapered tip facilitating easy introduction
- Good radiopacity for easy visualization
- The curvature of the S-shape dilator suits for the natural anatomy of male urethral

INTENDED USE

The Urethral Dilator is used for urethral dilation.

CLINICAL BENEFITS

In the operation of urethral stricture, urethral dilator is necessary to expand the urethral stricture to the normal urethral size, and the urethral dilator is placed into the urethral stricture from small size to large size in order to expand the urethral stricture to the normal urethral size. The clinical benefit associated with the device is define as direct, it is intended used urethral dilation. Prevention and treatment of inflammatory, traumatic and post-operative urethral stricture, and relief of urethral obstruction in child and adult.

INDICATIONS

Prevention and treatment of inflammatory, traumatic and post-operative urethral stricture, and relief of urethral obstruction.

PATIENT TARGET GROUP

Child and adult.

INTENDED USER

Must be operated by trained professionals.

CONTRAINDICATIONS

- Patients with acute urethritis, acute prostatitis or urethral purulent secretion.
- Patients with acute urethral injury, not to aggravate the damage, bleeding or shock or not to create a pseudo-urethra.
- Patients with suspected urethral tumor.
- Patients with urethral fever after dilatation of urethral stricture or urethral probe examination.
- Patients unable to receive implantation due to severe stricture caused by urinary tract atresia or other reasons.

ADVERSE REACTIONS

- Urethral perforation.
- Bleeding, inflammation or edema caused by urethral stimulation or damage.
- Damage to other adjacent organs.
- Pseudo-urethral formation.
- Infection.

DIRECTIONS FOR USE

- Choose the appropriate models of the urethral dilator according to the operation and patient situation.
- Check if the package is in good condition and confirm the component is intact. Check the model and size on the product label to make sure the proper model and size is used.
- Introduce the flexible tip of the urological guide wire into the urethral meatus, and gently manipulate it beyond the obstruction and into the bladder.

NOTE: The urological guide wire can also be placed through a cystoscope and maintained in position as the scope is removed.

- For standard type and standard S-curve type, lubricate the dilators with water-base lubricant. For hydrophilic type and hydrophilic S-curve type, soaking the dilators into sterile water to activate the coating.
- Pass the dilators over the urological guide wire while maintaining the urological guide wire's position, progressing from the smallest to the largest appropriate size.
- Remove the dilators and urological guide wire slowly once the dilation is completed.

WARNINGS

- Read all warnings and instructions before use. Improper use can result in serious or fatal illness or injury.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Do not use if it's damaged or irregularly shaped.
- Do not use after the expiry date.
- For single use only.
- Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the

patient.

PRECAUTIONS

- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.
- Use of this device should be restricted to use by or under the supervision of physicians trained.
- Performing dilation using dilators from small to large sizes in a gradual manner.
- Be slow and gentle when advancing and withdrawing the dilator.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the user's place.

STORAGE CONDITIONS

- Store in a cool and dry place, keep away from sunlight.
- Protect product from moisture and excessive heat.
- Store in manner preventing crushing.

APPENDIX A The Key Performance Characteristics of Urethral Dilator

The key performance characteristics: tensile strength

Test method: according to Annex C of ISO 11070:2014 A1:2018

Performance Requirement: the peak tensile force shall be $\geq 15\text{N}$.

MEANING OF SYMBOLS ON PACKAGE

Date of
manufacture

Use-by-date

Batch code

Do not re-use

Consult
instructions for
use

Do not use if
package is
damaged

**Sterilized
using ethylene
oxide**

Do not resterilize


**CE Marked
Product**

Medical Device

Manufacturer

**Authorized
representative
in the
European
Community**

Importer

Keep dry

**Keep away
from sunlight**

**Catalogue
number**

This way up

**Fragile, handle
with care**

**Temperature
limits:
-15°C~49°C**

**Single sterile
barrier system
(if applicable)**

**Single sterile
barrier system
with protective
packaging
inside
(if applicable)**

WELL LEAD MEDICAL CO., LTD.

**Address: C-4 Jinhu Industrial Estate, Hualong, 511434 Panyu, Guangzhou, PEOPLE'S
REPUBLIC OF CHINA**

E-mail: info@welllead.com.cn www.welllead.com.cn

Tel: +86 20 84758878

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg, GERMANY

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