PVC Laryngeal Mask Device

DESCRIPTION

The Qsafe Laryngeal Mask Device is an innovative, second generation, single use supraglottic airway management device.

The Qsafe Laryngeal Mask Device provides access to a functional separation of the respiratory and digestive tracts. The anatomically shaped airway tube is elliptical in cross section and ends distally at the laryngeal mask. The inflatable cuff is designed to conform to the contours of the hypopharynx, with the bowl and the mask facing the laryngeal opening—the Oropharyngeal Seal.

The Qsafe Laryngeal Mask Device also contains a drain tube which emerges as a separate port proximally and continues distally along the inside of the cuff bowl, passing through the distal end of the cuff to communicate distally with the upper oesophageal sphincter—the Oesophageal Seal. The drain tube may be used for the passage of a well lubricated gastric tube to the stomach, offering easy access for evacuation of gastric contents.


A built-in bite-block reduces the potential for tube damage and obstruction by patient biting. Attached to the mask is a cuff inflation line terminating in a pilot balloon and one-way check valve for mask inflation and deflation.

The Qsafe Laryngeal Mask Device is indicated for use in achieving and maintaining control of airway during routine and emergency anaesthetic procedures in fasted patients using either spontaneous or positive pressure ventilation.

The Qsafe Laryngeal Mask Device is made primarily of medical grade polyvinylchloride (PVC) and is supplied sterile for single use only. It consists of Cuff, Anatomically-shaped Airway Tube, Connector, Drain Tube, Cuff Inflation Line, Pilot Balloon, Valve, Bite-block.

The product can be used on adult patients.

Size	3#	4#	5#
Max. airway tube capacity, ml/cc	20	28	36
Max. Gastric tube size, Fr	14	14	14
Max. airway tube pressure drop, cmH2O	3	3	3
Distance from connector machine end to the ventilator opening, L	125	140	155

mm			
----	--	--	--

FEATURES

- Provides access to and functional separation of the respiratory and digestive tracts.
- Enables passive drainage or active management of digestive tract contents independent of ventilation.
- Designed to reduce gastric insufflation during ventilation and reduce the risk of aspiration.
- Bite-block helps to reduce potential dental damage.
- Elliptical and anatomically-shaped airway tube facilitates insertion success
- Soft & elongated cuff designed to support an effective seal
- Latex free

INTENDED USE

It is suitable to be used on patients under anesthesia or drug sedation and patients who need emergency access to artificial ventilation during first aid and resuscitation to achieve patency of the upper respiratory tract.

INDICATIONS

1. Sustained ventilation in the surgery under general anesthesia: Except the esophageal or gastrointestinal surgery with the risk for aspiration, the pulmonary or tracheal surgery with airway operation and the surgery for one-lung ventilation, as long as appropriate type and model are selected, laryngeal mask can be nearly used in any other surgeries. Laryngeal mask is particularly applicable to the following circumstances.

(1) Short body surface or limb surgery under general anesthesia: Ordinary laryngeal mask can be used to maintain autonomous-respiration general anesthesia or short-time IPPV general anesthesia. The second generation of SAD is recommended for direct IPPV general anesthesia.

(2) Head & neck surgery and ophthalmic surgery: Laryngeal mask is applicable to the head & neck minor surgery under general anesthesia, including ophthalmic surgery, E.N.T. surgery and plastic surgery. Flexible laryngeal mask can be used to reduce the effect on the surgical field. When intraocular surgery on the patient with elevated intraocular pressure, the insertion and

extraction of laryngeal mask have little effect on the intraocular pressure.

(3) Respiratory medicine and thoracic surgery: Insert the laryngeal mask under surface anesthesia plus sedation or general anesthesia, and retain spontaneous respiration by intravenous anesthesia or inhalation anesthesia for the following operations: ① fibrolaryngoscopy and fiberoptic bronchoscopy with laryngeal mask; ② endotracheal and supraclavicular tumor removal by Nd-YAG laser with laryngeal mask; ③ placement of trachea and bronchial dilator with laryngeal mask; ④ placement of fiberoptic bronchoscope with laryngeal mask in ICU, percutaneous dilatational tracheostomy under guidance of fiberoptic bronchoscope.

(4) Anesthesia with intra-operative wake-up test: Patients to be waken up sometimes during corrective surgery for scoliosis and resection of brain lesion in brain functional region are well-tolerated in the use of laryngeal mask, with shorter wake-up time, satisfactory in the cooperation, with little effect of position.

(5) Surgery in the lateral or prone position: Insertion of laryngeal mask is easier endotracheal intubation in the lateral or prone position. Prior to induction under general anesthesia, the patient can choose his/her own comfortable position. This is one of the advantages for the patients in the prone position when use of laryngeal mask. In addition, even if regurgitation, aspiration will not occur easily. The beginner is suggested to master the management of laryngeal mask in the prone position, and then to use the laryngeal mask when surgery in the lateral or prone position.

(6) Anesthesia outside the operating room: Some treatment or diagnosis operations for adults and children are performed outside the operating room, with little stimulation and mild pain, but the patients shall not move their positions, e.g. radiotherapy, angiography, interventional therapy, endoscopy and electrical cardioversion. Under intravenous induction, insert the laryngeal mask, retain the spontaneous respiration, and maintain intravenous anesthesia for the above operations.

2. Expected difficult airway. Those with expected difficulty in endotracheal intubation may directly use the laryngeal mask to maintain ventilation and anesthesia, and most of the patients with difficulty in endotracheal intubation may have the laryngeal mask inserted smoothly. When failure of endotracheal intubation, laryngeal mask may also be used as the ventilation tool or used for guiding the endotracheal intubation. However, such patients should guard against the risk of failure of intraoperative laryngeal mask ventilation. Therefore, the use of laryngeal mask shall be managed by an experienced anaesthetist.

3. Unexpected difficult airway. For the patients with induction under general anesthesia, when failure of endotracheal intubation, laryngeal mask may be used as the ventilation tool or used for guiding the endotracheal intubation.

4. Emergency airway. When difficulty in mask ventilation, particularly when "neither endotracheal intubation nor mask ventilation", establish effective

Instruction for Use

ventilation and oxygenation as soon as possible. Laryngeal mask is often used as the preferred tool. At this time, laryngeal mask may be used for ventilation, or endotracheal intubation may be performed through the laryngeal mask. Other supraglottic ventilation tools (e.g. venturi) may also be used for emergency airway.

5. Tracheal extubation. When insertion of laryngeal mask before or after extubation, or after tracheal extubation, laryngeal mask may be used for ventilation or used for guiding the endotracheal intubation when necessary. Replacement of tracheal tube with laryngeal mask may establish a physiologically safe and non-stimulated airway, which can prevent the secretion and blood from oral cavity. This technology may be used for wake-up extubation or extubation under deep anesthesia, particularly applicable to the patients with cardiovascular system stimulation caused by tracheal tube that might affect the effect of surgical repair and the patients with coronary heart disease or other heart diseases, and possibly more beneficial to the patients with smoking, asthma or other airway hypersensitivity.

6. First-aid and resuscitation in the prehospital environment, emergency department, ICU and other departments.

CLINICAL BENEFITS

As life support to establish artificial airway, ensure circulation, ventilation and oxygen supply, and quickly recover and maintain effective respiration of patients.

PATIENT TARGET GROUP

Adult patients.

INTENDED USER

Must be operated by trained professionals.

CONTRAINDICATIONS

Based on the current development trend, the indications and contraindications of laryngeal mask are constantly changing. Generally, the contraindications are mostly relative contraindications, which shall be taken into account according to the actual condition of the patient, the criticality of airway treatment and the special requirements for surgical operation.

1. Patients at high risk of gastric content reflux and respiratory aspiration, including non-fasting, full stomach, morbid obesity, frequent gastroesophageal reflux, intestinal obstruction and esophageal hiatal hernia.
2. Patients with tumor, injury, infection or other pathological changes in the laryngeal part of pharynx, e.g. supraglottic or hypopharyngeal injury, tumor, severe tonsil hypertrophy, significant laryngeal or tracheal deviation, softening and peripheral compression.
3. Patients with decreased pulmonary compliance, particularly those with expected intraoperative airway pressure greater than 30cmH₂O.

4. Those with difficulty in insertion of laryngeal mask due to small mouth opening.
5. Those with expected long surgery time, although it was once reported that laryngeal mask could be used safely for more than 8h.

ADVERSE REACTIONS

- Aspiration pneumonia
- Airway obstruction
- Laryngospasm
- Sore throat
- Dryness of the throat and/or mucosa
- Laryngeal nerve injury and vocal cord paralysis
- A diminished sense of taste
- Side effects due to improper placement vary based on the nature of the placement.


DIRECTIONS FOR USE

- Choose proper Qsafe Laryngeal Mask Device size according to actual situation of patient.
- Check the model and size on the product label to make sure the proper model and size is used.
- Examine the Qsafe Laryngeal Mask Device before use to make sure that the device and the packaging is not damaged or without leaking.
- Open the package and take out the device. Use syringe to deflate the cuff.
- Lubricate the posterior surface of the mask and airway tube just prior to insertion.
- The operator gently holds the jaw to assist the patient to open the mouth. Insert the mask along the back of the tongue.
- Press the tip of the mask against the hard palate. Press the cuff further into the mouth, maintaining pressure against the palate. Swing the device inward with a circular motion, pressing against the contours of the hard and the soft palate. Advance the device into the hypopharynx until resistance is felt.
- Correct placement should produce a leak-free seal against the glottis with the mask tip at the upper oesophageal sphincter. The integral bite block should lie between the teeth.
- The device placement should be confirmed again if the patient's body position is changed.

Table 1: Selection Guide

Size	Applicable Patient Weight kg	Recommended Maximum Inflation Volume ml	Largest size of placed gastric tube Fr	Maximum outer diameter of device passing the Laryngeal Mask mm	Optimum Intra-Cuff Pressure cmH ₂ O
3#	30-50	30	14	5.3	60
4#	50-70	40	14	6	
5#	70-100	40	14	6	

1 esophagus, 2 trachea, 3 vocal cord/folds, 4 thyroid cartilage, 5 laryngeal inlet, 6 epiglottis, 7 hyoid bone, 8 tongue, 9 incisors

WARNINGS

- Read all warnings and instructions before use. Improper use can result in serious or fatal illness or injury.
- This device can only be used by trained professionals.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Do not use if it's damaged or irregularly shaped.
- Do not use after the expiry date.
- Not compatible with MRI.

PRECAUTIONS

- When applying lubricant, prevent blocking the airway
- Suggest maximum cuff pressure 60 cm H₂O
- Diffusion of nitrous oxide (N₂O) mixture, oxygen or air from a surrounding gas mixture may either increase or decrease cuff volume and pressure.
- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.
- Be care for the risk of inflammable if use any laser device.

Instruction for Use

- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the user's place.

STORAGE CONDITIONS

- Store in a cool and dry place, keep away from sunlight.
- Protect product from moisture and excessive heat.
- Avoid prolonged exposure to ultraviolet, sunlight and fluorescent light.
- Store in manner preventing crushing.

The Key Performance Characteristics of the Qsafe Laryngeal Mask Device

Key Performance: Kink Resistance

The pressure drop of the airway complies with the requirements in the table 1.

Table 1 Maximum Pressure Drop

Size	3#	4#	5#
Maximum Pressure Drop, cmH2O	3	3	3

MEANING OF SYMBOLS ON PACKAGE

	Date of manufacture		Use-by-date		Batch code
	Do not re-use		Consult instructions for use		Do not use if package is damaged
	Sterilized using ethylene oxide		Do not resterilize		Single sterile barrier system
	Medical Device		Manufacturer		Authorized representative in the European Community
	CE Marked Product		Keep dry		Keep away from sunlight

wellead

Instruction for Use

Importer

This way up

Fragile, handle with care

Not contains or presence of phthalate.

Catalogue number

Temperature limits: -15°C~49°C

WELL LEAD MEDICAL CO., LTD.

Address: C-4 Jinhua Industrial Estate, Hualong, 511434 Panyu, Guangzhou, PEOPLE'S REPUBLIC OF CHINA

E-mail: info@welllead.com.cn www.welllead.com.cn

Tel: +86 20 84758878

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg, GERMANY

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Shanghai International Holding Corp. GmbH (Europe) Address: Eiffestrasse 80, 20537 Hamburg, GERMANY